

eCOA Solutions

Interactive Wireframe Guide

A comprehensive guide to the UI/UX design and screen flows

Version 1.0 | December 2025

1. Overview

This document describes the interactive wireframe prototype for eCOA Solutions. The prototype demonstrates the core user flows for sponsor administrators to create and manage clinical trial forms.

Key Features:

- Modern, dark-themed UI with cyan/violet accent colors
- Responsive design optimized for desktop and tablet
- Intuitive drag-and-drop form builder
- Real-time properties panel for element customization
- Multi-screen navigation with smooth transitions

2. Screen Descriptions

2.1 Login Screen

The entry point for sponsor administrators to access the platform.

Element	Description
Logo	eCOA Solutions branding with gradient icon
Username Field	Text input with placeholder and validation
Password Field	Secure input with masked characters
Sign In Button	Gradient button (cyan → violet) with hover effect
Background	Dark gradient with subtle pattern overlay

2.2 Dashboard Screen

The main hub displaying all studies for the logged-in sponsor.

Element	Description
Sidebar	Collapsible navigation with Dashboard, Studies, Patients, Reports, Settings
Header	'My Studies' title with 'Create New' button
Search Bar	Filter studies by name
Study Cards	Show name, status badge (Draft/UAT/Production), last modified date
Status Badges	Color-coded: Green (Production), Amber (UAT), Gray (Draft)

2.3 Create Study Modal

Modal dialog presenting study creation options.

Option	Description
Use Templates	Opens template selection with Clone Study, Licensed Quest, Company Templates, C
Create from Scratch	Opens blank form builder canvas for custom form creation

2.4 Template Selection Screen

Browse and select from available form templates.

Tab	Content
Clone Study	Search existing studies to duplicate
Licensed Quest	PHQ-9, SF-36, PROMIS validated instruments
Company Templates	Organization-specific templates (e.g., 'Cipla Templates')
Our Templates	Custom saved templates

2.5 Form Builder Screen

The core form creation interface with drag-and-drop functionality.

Panel	Description
Tools Panel (Left)	Draggable elements: Radio, Checkbox, Scale 1-3, Text, Add Logo
Canvas (Center)	Drop zone for form elements, shows 'Blank Canvas' when empty
Properties Panel (Right)	Font, Size, Color, Alignment settings for selected element
Top Bar	Form name input, Preview button, Finalize button
Bottom Bar	Back/Next navigation, page indicator

2.6 Form Properties Screen

Configure form settings before deployment.

Setting	Options
Completed By	Patient, Caregiver, Investigator (multi-select)
Notification Settings	Frequency (Once/Twice daily), Time picker
Schedule & Triggers	Daily, Weekly, Monthly, Custom
Include in Reports	Toggle switch (Yes/No)
Actions	Save as Draft, Deploy to UAT buttons

3. User Flows

3.1 Create Study from Template

1. Login → Dashboard
2. Click 'Create New' button
3. Select 'Use Templates' option
4. Choose template category (Clone/Licensed/Company/Custom)
5. Search and select specific template
6. Click 'Clone' to open in Form Builder
7. Customize form as needed
8. Click 'Finalize' to configure properties
9. Click 'Deploy to UAT' to launch

3.2 Create Study from Scratch

1. Login → Dashboard
2. Click 'Create New' button
3. Select 'Create from Scratch' option
4. Form Builder opens with blank canvas
5. Drag elements from Tools panel to canvas
6. Click element to edit properties
7. Add multiple pages as needed
8. Click 'Finalize' to configure settings
9. Click 'Deploy to UAT' to launch

4. Technical Implementation

The wireframe is built as an interactive React component using:

- React 18 with functional components and hooks
- Tailwind CSS for styling (dark theme)
- Lucide React for iconography
- State management via useState hooks
- Component-based architecture for reusability

To run the prototype:

1. Copy eCOA_Interactive_Wireframe.jsx to a React project
2. Install dependencies: `npm install lucide-react`
3. Import and render the component
4. Or open directly in Claude.ai to see the interactive preview